



**DALMIA CEMENT (B) LIMITED – CEMENT PLANT
ARIYALUR
INSTRUMENTATION MANUAL**

Issue No. 02	Rev. No: 00	Effective Date: 01.11.2019	INST/SOP/016
Issued By: M.R		Approved By: HOD - INSTRUMENTATION	
SOP for Maintenance of ROBOT Lab Equipments.			

Scope : This SOP is applicable for the maintenance of robot lab equipments.

Responsibility : Instrument Technician.

Accountability : Instrument -Section Engineer.

PPE:

1. Safety helmet,
2. Safety shoe,
3. Mask.
4. Gloves.
5. Safety googles.

Tools:

- Multimeter
- Tester
- Insulated Screw Driver sets
- Hand Blower
- Ring spanner sets

Hazards:

Risks associated

Electrocution

Mitigating Measures

Use Insulated hand glove and Tools

Training needs:

1. People should have knowledge to work in Robot lab.
2. Emergency procedure



**DALMIA CEMENT (B) LIMITED – CEMENT PLANT
ARIYALUR
INSTRUMENTATION MANUAL**

Issue No. 02	Rev. No: 00	Effective Date: 01.11.2019	INST/SOP/016
Issued By: M.R		Approved By: HOD - INSTRUMENTATION	
SOP for Maintenance of ROBOT Lab Equipments.			

Procedure:

1. Use Proper PPE & tools.
2. Take work permit & Line clearance for Robot lab particular equipments (Receiving station, Mill & press 1&2, Bottle magazine, Manual drawer, Cup cleaning, Ring cleaning, Robot, Sending stations).
3. Ensure that Robot in rest position and manual mode. If we are working in Robot movement area.
4. If robot is not in rest position after reaching the home position put in manual mode by using Auto/Manual Key.
5. Close the incoming compressor Air for particular Equipment and till wait for accumulated air vent from equipment valves & cylinders by removing the outlet port.
6. Switch off the Isolator / MCB.
7. Isolate the control supply by switching off the MCB.
8. Ensure non availability of incoming supply at incoming terminal by using multimeter.
9. Ensure the Earthing connection of panels.
10. Clean the Panel with hand blower / Vacuum cleaner.
11. Check the tightness of the connections by using proper insulated tools.
12. Check the tightness of foundation bolts.
13. Clean /Check the solenoid valve silencers.
14. After the panel maintenance ensures the panel Screws and door locks are properly locked.
15. Switch on the Isolators /MCB and control supply MCB.
16. Open the compressor air line and arrest if there is any air leakage.
17. Reset the equipment by pressing the reset equipment.
18. Ensure the healthiness of the equipment.
19. Ensure the robot home position and Auto mode. If we are worked in Robot movement area.
20. Give the clearance for operation



**DALMIA CEMENT (B) LIMITED – CEMENT PLANT
ARIYALUR
INSTRUMENTATION MANUAL**

Issue No. 02	Rev. No: 00	Effective Date: 01.11.2019	INST/SOP/016
Issued By: M.R		Approved By: HOD - INSTRUMENTATION	
SOP for Maintenance of ROBOT Lab Equipments.			

Emergency / Emergency Shut OFF:

1. If any unconsciousness is there, give First aid and inform to the Safety department or Call Emergency number 3108/9865152222/9750963761.

Records/Annexure:

1. Refer Line clearance certificates.
2. JSA as enclosed below.

JOB SAFETY ANALYSIS: (JSA)

Job Safety Analysis	Job: Robot lab	Date: 20.9.2017	Analysis by: Section Head	Reviewed by: HOD
Title of employee doing job: INST - TECHNICIAN	Supervisor: Section Engineer	Department: Instrumentation	Section: Robo lab	Approved by:
Req'd/recommended PPE: Safety helmet, Safety shoe, Mask, Goggles & Gloves.				
Sequence of Basic Job Steps	Potential Hazards	Recommended Safe Job Procedure	What Could Go Wrong	Corrective Action
1.TAKE THE LINE CLEARANCE & SWITCH OFF THE ISOLATOR	LIVE SUPPLY 3 PHASE IS EXPOSED AT THE PANEL SIDE	ENSURE THE ISOLATION OF POWER SUPPLY AT THE INCOMING SIDE & CHECK THE INCOMING SUPPLY WITH MULTIMETER.	RESULT IN ELECTROCUT ED& FLASHING.	CHECK THE INCOMING SUPPLY WITH MULTIMETER.



**DALMIA CEMENT (B) LIMITED – CEMENT PLANT
ARIYALUR
INSTRUMENTATION MANUAL**

Issue No. 02	Rev. No: 00	Effective Date: 01.11.2019	INST/SOP/016
Issued By: M.R		Approved By: HOD - INSTRUMENTATION	
SOP for Maintenance of ROBOT Lab Equipments.			

2.ENSURE THE ROBOT IN REST POSITION	ROBOT MAY HIT THE PERSON	ENSURE THE ROBOT REST POSITION AND ROBOT KEY IN MANUAL MODE.	RESULT IN PERSON INJURY	CHECK THE ROBOT POSITION AND MANUAL MODE IN PENDANT
3.CLOSE THE COMPRESSOR AIR LINE	POSSIBLE FINE DUST AND RUST PARTICLES COME FROM HOSE LINE.	WEAR GOOGLES, HAND GLOVES & NOSE MASK	COUGHING, EYE IRRITATION MAY OCCUR	USE GOOGLES, DUST MASK AND HAND GLOVES
3.CLEANING OF PANEL	POSSIBLE FINE DUST COME OUT	WEAR GOOGLES, HAND GLOVES & NOSE MASK	COUGHING, EYE IRRITATION MAY OCCUR	USE GOOGLES, DUST MASK AND HAND GLOVES
4.CHECK THE CONTROL CABLE TIGHTNESS	CHANCES FOR TERMIANL BLOCK DAMAGE	USE PROPER TOOLS BASED ON THE TYPES OF TERMINAL BLOCKS	TERMINALS GET DAMAGE AND NUISENCE OPERATION DURING PANEL STARTUP	SELECT THE PROPER TOOLS
5.CHECK THE POWER CABLE TIGHTNESS	CHANCES FOR TERMINAL BLOCK DAMAGE	USE PROPER TOOLS BASED ON THE TYPES OF TERMINAL BLOCKS	TERMINALS GET DAMAGE AND NUISENCE OPERATION DURING PANEL STARTUP	SELECT THE PROPER TOOLS



**DALMIA CEMENT (B) LIMITED – CEMENT PLANT
ARIYALUR
INSTRUMENTATION MANUAL**

Issue No. 02	Rev. No: 00	Effective Date: 01.11.2019	INST/SOP/016
Issued By: M.R		Approved By: HOD - INSTRUMENTATION	
SOP for Maintenance of ROBOT Lab Equipments.			

6.CHECK THE TIGHTNESS OF BOLTS & NUTS	CHANCES FOR BOLTS & NUTS HEAD MAY BE DAMAGE	USE PROPER TOOLS BASED ON THE TYPES OF BOLTS & NUTS	BOLTS & NUTS GET DAMAGE DURING OPEARTION EQUIPMENTS MAY BE CREATING NUISENCE NOISE.	SELECT THE PROPER TOOLS
7.CHECK THE TIGHTNESS OF HOSE CLIPS	CHANCES FOR BOLTS & NUTS HEAD MAY BE DAMAGE	USE PROPER TOOLS BASED ON THE TYPES OF BOLTS & NUTS	BOLTS & NUTS GET DAMAGE DURING OPEARTION EQUIPMENTS MAY BE CREATING NUISENCE NOISE OR AIR LEAKAGE.	SELECT THE PROPER TOOLS